


Facial Implants




sientra®

Facial Implants

The Sientra Facial Implants are manufactured from medical grade silicone, for the purpose of facial reconstructive or cosmetic interventions. A large variety of models and formats is available to meet surgeon needs. Supplied sterile.


Nasal "L" Implant

Smooth Surface



SHORE	DIMENSIONS				REFERENCE
	A cm	B cm	C cm	D cm	
25/50	1.90	4.10	2.80	0.40	71-060

It has two different silicone hardnesses, being more rigid in the dorsum. It is radiopaque and may be trimmed as needed.



SHORE	DIMENSIONS				REFERENCE
	A cm	B cm	C cm	D cm	
10/40	1.80	4.00	2.90	0.45	72-060


It has two different silicone hardnesses, being more rigid in the dorsum. It is radiopaque and may be trimmed as needed.



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Anatomical Solid Chin Implant

Smooth Surface




SIZE	DIMENSIONS			REFERENCE
	A cm	B cm	C cm	
1= Small	6.10	1.40	0.60	15-060
2= Medium	6.40	1.50	0.70	15-065
3= Large	7.90	1.50	0.80	15-080

Anatomical chin implants have a reinforced patch for fixation stitches.

Round Solid Chin Implant

Smooth Surface




SIZE	DIMENSIONS			REFERENCE
	A cm	B cm	C cm	
1= Small	3.50	1.60	0.70	10-035
2= Medium/ Small	4.50	1.70	0.80	10-045
3= Medium/ Large	4.80	1.70	0.80	10-050
4= Large	5.30	1.80	0.80	10-055
5= Trimmable	7.50	1.90	0.80	10-075

Round chin implants have holes to encourage tissue growth.

Anatomical Malar Implant

Smooth Surface



SIZE	DIMENSIONS			REFERENCE
	A cm	B cm	C cm	
1= Small	2.00	4.00	0.42	25-40D
	2.00	4.00	0.42	25-40E
2= Medium	2.30	4.30	0.42	25-45D
	2.30	4.30	0.42	25-45E
3= Large	2.50	5.00	0.45	25-50D
	2.50	5.00	0.45	25-50E

Anatomical malar implants have a lateral extension and holes to encourage tissue growth.




Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Nostril Retainer


The Sientra Nostril Retainer was designed for prophylaxis of stenosis and cicatricial contractures. It is made of molded silicone elastomer, allowing its length to be adjusted to each patient's particular need. It has a functional, anatomic format, two reinforced silicone strips for fixation as well as a white dot to indicate proper positioning. Designed to maintain the septum in vertical position to allow for nasal breathing immediately following the operation without blocking the flow of air. Supplied sterile in individual packages.

INDICATIONS:

- After aesthetic and reconstructive rhinoplasty
- After rhinoseptoplasty
- After rhinoseptocheiloplasty and/or primary and secondary rhinocheiloplasty (cleft-lip nose patients)



CORRECT POSITION



SIZE	DIMENSIONS			REFERENCE
	A cm	B cm	C cm	
1	1.80	2.30	0.70	1171-018
2	1.90	2.40	0.80	1171-019
3	2.00	2.50	0.90	1171-020
4	2.10	2.60	1.00	1171-021
5	2.20	2.70	1.05	1171-022
6	2.30	2.80	1.10	1171-023
7	2.40	2.90	1.20	1171-024
8	2.50	3.00	1.30	1171-025
9	2.60	3.10	1.40	1171-026
10	2.70	3.20	1.50	1171-027
11	2.80	3.30	1.60	1171-028
12	2.90	3.40	1.70	1171-029
13	3.00	3.50	1.80	1171-030



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.