

# Gluteal Implants




sientra.

# Gluteal Implants

Sientra gluteal implants are composed of a silicone elastomer envelope filled with silicone elastomer. They are used in reconstructive and cosmetic surgery, and provided in three distinct shapes. The Oval Shaped implants feature an orientation line for positioning.

## Oval Gluteal Implant


Textured Surface



DIMENSIONS				REFERENCE
A cm	B cm	C cm	VOL cc	
15.3	10.3	2.3	225 cc	20922-150
15.3	10.3	2.5	260 cc	20922-155
15.6	10.9	2.7	300 cc	20922-160
16.0	11.5	2.9	330 cc	20922-165
16.3	11.9	3.2	385 cc	20922-168
16.6	12.8	3.3	420 cc	20922-170

## Oval Gluteal Implant

Textured Surface




DIMENSIONS				REFERENCE
A cm	B cm	C cm	VOL cc	
13.6	10.0	2.0	145 cc	20942-135
14.5	10.4	2.9	233 cc	20942-145
16.0	11.6	3.0	275 cc	20942-160
17.1	12.4	3.2	343 cc	20942-170
17.3	12.8	3.5	385 cc	20942-175
17.8	13.2	3.9	435 cc	20942-180
18.4	13.8	4.1	540 cc	20942-185

## Round Gluteal Implant

Smooth Surface



DIMENSIONS				REFERENCE
A cm	B cm	C cm	VOL cc	
11.6	4.0	2.6	250 cc	932-115
12.1	4.3	2.9	290 cc	932-120
12.5	4.5	2.9	340 cc	932-125
13.3	5.5	4.0	440 cc	932-130



### GLUTEAL IMPLANT

- 1) IMPLANT
- 2) GLUTEAL MEDIUS
- 3) GLUTEAL MAXIMUS
- 4) SCIATIC NERVE



**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.