

QUICK REFERENCE GUIDE

teardrop shape

HSC+ textured breast implants & saline sizers

see your skills in
sientra


Doctor _____ Surgery Date _____ Patient Name _____

Phone Number _____ Billing Zip Code _____ Email _____

Office Contact _____ Surgery Location _____


Comments _____

Classic Base Moderate Projection




WTH cm	HT cm	PROJ cm	VOL cc	HSC+		SIZER	
				Catalog #	QTY	Catalog #	QTY
9.5	10.5	3.6	180cc	CBASE-180MP		SZ20676-180E	
10.0	11.0	3.8	210cc	CBASE-210MP		SZ20676-210E	
10.5	11.5	4.0	240cc	CBASE-240MP		SZ20676-240E	
10.8	11.8	4.1	260cc	CBASE-260MP			
11.0	12.0	4.2	275cc	CBASE-275MP		SZ20676-275E	
11.3	12.3	4.3	295cc	CBASE-295MP			
11.5	12.5	4.4	310cc	CBASE-310MP		SZ20676-310E	
11.8	12.8	4.5	335cc	CBASE-335MP			
12.0	13.0	4.6	350cc	CBASE-350MP		SZ20676-350E	
12.3	13.3	4.7	375cc	CBASE-375MP			
12.5	13.5	4.8	400cc	CBASE-400MP		SZ20676-400E	
12.8	13.8	4.9	425cc	CBASE-425MP			
13.0	14.0	5.0	450cc	CBASE-450MP		SZ20676-450E	
13.3	14.3	5.1	475cc	CBASE-475MP			
13.5	14.5	5.2	500cc	CBASE-500MP			
13.8	14.8	5.3	525cc	CBASE-525MP			
14.0	15.0	5.4	550cc	CBASE-550MP			
14.3	15.3	5.5	575cc	CBASE-575MP			
14.5	15.5	5.6	600cc	CBASE-600MP			
14.7	15.7	5.7	625cc	CBASE-625MP			
14.9	15.9	5.7	650cc	CBASE-650MP			
15.0	16.0	5.8	700cc	CBASE-700MP			

Classic Base High Projection




WTH cm	HT cm	PROJ cm	VOL cc	HSC+		SIZER	
				Catalog #	QTY	Catalog #	QTY
9.0	10.0	4.2	190cc	CBASE-190HP			
9.5	10.5	4.4	215cc	CBASE-215HP			
10.0	11.0	4.6	250cc	CBASE-250HP			
10.5	11.5	4.8	285cc	CBASE-285HP			
10.8	11.8	4.9	300cc	CBASE-300HP			
11.0	12.0	5.0	320cc	CBASE-320HP			
11.3	12.3	5.1	345cc	CBASE-345HP			
11.5	12.5	5.2	365cc	CBASE-365HP			
11.8	12.8	5.3	390cc	CBASE-390HP			
12.0	13.0	5.4	410cc	CBASE-410HP			
12.3	13.3	5.5	440cc	CBASE-440HP			
12.5	13.5	5.6	460cc	CBASE-460HP			
12.8	13.8	5.7	485cc	CBASE-485HP			
13.0	14.0	5.8	510cc	CBASE-510HP			
13.3	14.3	5.9	545cc	CBASE-545HP			
13.5	14.5	6.0	570cc	CBASE-570HP			
13.8	14.8	6.1	600cc	CBASE-600HP			
14.0	15.0	6.2	635cc	CBASE-635HP			

Round Base Moderate Projection



WTH cm	PROJ cm	VOL cc	HSC+		SIZER	
			Catalog #	QTY	Catalog #	QTY
9.2	4.0	160cc	RBASE-160MP			
9.5	4.1	170cc	RBASE-170MP			
9.9	4.2	190cc	RBASE-190MP			
10.1	4.3	205cc	RBASE-205MP			
10.4	4.4	220cc	RBASE-220MP			
10.7	4.5	245cc	RBASE-245MP			
11.0	4.6	260cc	RBASE-260MP			
11.3	4.7	275cc	RBASE-275MP			
11.6	4.8	300cc	RBASE-300MP			
11.9	4.9	320cc	RBASE-320MP			
12.2	5.0	340cc	RBASE-340MP			
12.5	5.1	360cc	RBASE-360MP			
12.8	5.2	390cc	RBASE-390MP			
13.1	5.3	415cc	RBASE-415MP			
13.4	5.4	445cc	RBASE-445MP			
13.7	5.5	470cc	RBASE-470MP			
14.0	5.6	500cc	RBASE-500MP			
14.3	5.7	540cc	RBASE-540MP			
14.6	5.8	595cc	RBASE-595MP			
14.9	5.9	625cc	RBASE-625MP			
15.2	6.0	650cc	RBASE-650MP			
15.5	6.1	700cc	RBASE-700MP			

Round Base High Projection



WTH cm	PROJ cm	VOL cc	HSC+		SIZER	
			Catalog #	QTY	Catalog #	QTY
9.8	4.9	200cc	RBASE-200HP			
9.9	5.0	220cc	RBASE-220HP		SZ20646-220RB	
10.1	5.1	235cc	RBASE-235HP		SZ20646-235RB	
10.4	5.2	255cc	RBASE-255HP		SZ20646-255RB	
10.7	5.3	270cc	RBASE-270HP		SZ20646-270RB	
11.0	5.4	295cc	RBASE-295HP			
11.3	5.5	320cc	RBASE-320HP		SZ20646-320RB	
11.6	5.6	345cc	RBASE-345HP			
11.9	5.7	370cc	RBASE-370HP		SZ20646-370RB	
12.2	5.8	395cc	RBASE-395HP			
12.5	5.8	425cc	RBASE-425HP		SZ20646-425RB	
12.8	6.0	455cc	RBASE-455HP			
13.1	6.1	485cc	RBASE-485HP			
13.4	6.2	515cc	RBASE-515HP			

To Order Call 1.888.708.0808